

EVIS LUCERA DUODENOVideoscope

TJF-260V

Breakthrough V-Groove Design Ensures Fast,
Easy ERCP Device Replacement and Secure Device Manipulation



Reliable Device Replacement, Secure Guidewire Positioning, and Easier Handling

Equipped with OLYMPUS' innovative V-Groove forceps elevator, the revolutionary TJF-260V makes it possible to perform ERCP device manipulation and replacement with confidence and efficiency. The unique V-Groove on the forceps elevator securely locks the guidewire into position.



V-SYSTEM™

V-System™ device replacement procedure

Confirm the position of the V-Marking on the Endo-Therapy accessory.



When the V-Marking is completely visible above the instrument channel port, raise the forceps elevator to lock the guidewire.



Be sure that the guidewire is locked in the V-Groove.




Completely remove the device.



V-Groove forceps elevator

The innovative design of the forceps elevator on the TJF-260V allows it to be raised at a sharper angle than a conventional scope elevator, while its unique V-shaped guidewire locking groove combines with the dedicated guidewire's reactive force to securely lock the guidewire in position.



OLYMPUS TJF TYPE 260V



Main Features

- 4.2 mm instrument channel for expanded treatment options in the bile and pancreatic ducts.
- Innovative V-Groove forceps elevator offers reliable duodenoscopic treatment capability in combination with dedicated V-System ERCP devices.
- Insertion tube measuring 11.3 mm across with distal end diameter of 13.5 mm.
- Four-way angulation (120° Up, 90° Down, 110° Right, and 90° Left) facilitates approach to the papilla of Vater.
- Ergonomically designed grip enhances torque performance and scope maneuverability while easy-to-access control knobs improve operability.
- Detachable distal cover for easier cleaning of the scope tip.
- Fully compatible with the CV-260SL, CV-260 and CV-240.
- Scope ID function retains individual scope information in the memory chip and displays it on the monitor. Also stores settings such as auto white balance to facilitate endoscopy suite management.



Specifications

Optical System	Field of view	100°
	Direction of view	Backward oblique viewing 15°
	Depth of field	5 to 60 mm
Distal End	Outer diameter	13.5 mm
Insertion Tube	Outer diameter	11.3 mm
Bending Section	Angulation range	Up 120°, Down 90°, Right 110°, Left 90°
Working Length		1240 mm
Total Length		1550 mm
Instrument Channel	Inner diameter	4.2 mm
	Minimum visible distance	10 mm from distal end
	Endotherapy accessory entrance/exit position in field of view	
Laser Compatibility		NO
High-Frequency Compatibility		YES
Scope ID Function		YES



Distal end 13.5mm

Channel 4.2mm

CV Compatibility

CV-260SL
 CV-260
 CV-240



Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.