



Surgical Instruments

From the design to maintenance and care
and the function test...

Everything you need to know about our
instruments!



Table of Contents


Preface	4
KLS Martin Quality Instruments	6-7
General Information	
about the entire instrument range	8
Basic Principles	
of the maintenance, care and preparation of surgical instruments	9-12
Cutting	13
Scissors	14
Wire cutting pliers	15
Punches	16
Bone rongeur forceps	17
Osteotomes, chisels, gouges	18
Scalpels	19
Holding/Gripping	21
Clamps	22
Needle holders	23
Forceps	24
Pliers	25
Suctioning/Rinsing	27
Suction cannula	28
Trocar	29
Cannula	30
Syringe	31
Spreading	33
Wound spreaders	34
Abdominal retractors	35
Wound retractors	36
Rib spreader	37
Speculum	38

Preface

In order to maintain the value of high quality instruments, maintenance, care and preparation play a very important and central role.


This brochure aims to clarify basic questions about these points and give an overview of the diversity of high quality instruments provided by the KLS Martin Group and give the intended use of each of the product features.

OR



Storage


Disposal



Sterile Supplies Cycle

STERILE !

Sterilization



Ultrasound




Machine cleaning




Final rinsing with demineralized water




Manual cleaning



Care



1 Max. in the first lock




Functional test


Cleaning




All new instruments have to be washed at least three times before they can be included in the instrument cycle.




Single use



Defective instruments



Rusting instruments



Instruments for other sterilization procedures

Care and Functional Tests

KLS Martin Quality Instruments

Medical devices are subject to a set of uniform laws throughout the European Union. In this context, the EC Directive concerning medical devices places rigid quality demands on surgical instruments as well. Any manufacturer must heed and meet these requirements, as non-complying products are barred from being placed on the European single market. Similar regulations are in place in the United States and Japan. Other countries have started preparing and implementing their corresponding regulations.

Introduction

As an international supplier of medical devices, Gebrüder Martin GmbH & Co. KG has always given top priority to quality issues with regard to its surgical instruments. In fact, quality was a priority in our company long before legal regulations were put in place. Furthermore, since we have always gone a good deal beyond legal requirements, it is not surprising that competent KLS Martin employees have contributed greatly to the preparation and continuous improvement of both German and international standards governing the function, shape and quality of the various different instruments.

Development

It goes without saying that before an instrument can be manufactured, it must be designed first. The KLS Martin product divisions are, so to speak, the mediators and interpreters between the market and the design engineers. For any new instrument, they identify and define the features and properties desired by the user. In many cases, this involves bringing the development department in contact with leading surgeons, physicians or researchers. This collaboration is frequently reflected in the name given to the new product.

Once the requirements have been established, the design engineers can start developing the new instrument. Where new product series are concerned, we usually make use of our good relationships with leading research and testing institutes. All development efforts, including the tests performed and insights gained, are duly documented in the course of the design and development process. Upon completion of the development phase, the new product is validated against both the requirements defined by the product division and those defined by the law. A declaration of conformity is issued only for products that successfully passed these tests. The design is then released for production.

Materials

The materials used for the manufacture of surgical instruments are nationally and internationally standardized. As most instruments require a high mechanical strength for proper functioning, we are using hardenable chromium steels with a low to medium carbon content. In fact, a chromium content of at least 12.5% is required in order to guarantee sufficient corrosion resistance. As nickel-chromium steels, which have a distinctly higher corrosion resistance, are not heat-treatable, they can only be used for trays and specific large-surface instruments (see Annex "Current reference standards for KLS Martin products").

Raw Materials

The first step in manufacturing surgical instruments is forging a blank. There are only few certified specialized smithies that are capable of manufacturing such blanks in line with KLS Martin's specifications regarding the material, shape and dimensions.

Only stainless steel of European origin is used as a source material for manufacturing the blanks. Each steel batch is delivered with a recognized test certificate, and the blanks produced from this source material are likewise duly documented. The quality systems of these contractors have been well aligned with KLS Martin's own quality system. This guarantees that any product can be easily traced to the extent required by the law.

Manufacture

Based on the technical documentation established, the product then enters the manufacturing stage. All manufacturing specifications are subject to quality management. Any changes made on the product can therefore be easily tracked and pinned down at any time.

The manufacturing process is a batch process because the enormous variety of instruments required by our customers does not allow for serial production. This means that almost all instruments are, to a large degree, hand-made. Intermediate tests and inspections are clearly defined steps in the process and as such included in the manufacturing documentation.

Heat Treatment

From a functional and reprocessing perspective, the most important step in the manufacturing process is the heat treatment. The heat treatment gives those instruments that are made of hardenable chromium steels the required strength, tenacity and corrosion resistance. In contrast, instruments made of nickel-chromium steels are not hardenable, which means that the use of such steels is limited to special instruments.

In a first step, called "hardening" or "quenching", the instruments are heated to a temperature of more than 1,000°C (1,832°F). At this very high thermal level, the chromium-carbon compounds previously encapsulated in the material dissolve completely, distributing uniformly in the steel. To preserve this ideal structure, the steel is "frozen" by cooling it down fast, which results in a needle-type structure that gives the instrument the strength required for fulfilling its function. At the same time, the uniform distribution of the chromium content enhances the corrosion resistance to such a degree that the instruments are fully corrosion-resistant as long as the user respects the conditions specified for instrument processing.

In a second step, called "tempering", the instruments are kept at a temperature of about 250°C (482°F) for several hours. This treatment reduces the tensions present in the instrument. As a result, the tempered instruments are significantly more elastic and less prone to fractures.



Proper heat treatment is directly reflected in the measured hardness values (see Annex “Current reference standards for KLS Martin products“). Lower (substandard) hardness values indicate insufficient heat treatment and, therefore, a lack of functionality and corrosion resistance.

Function and Finish

Following the heat-treatment stage, specially trained instrument makers give the instruments their final shape, function and finish.

Notably, Tuttlingen is the only place in the world with a system of technical colleges and training institutions specially adapted to the needs associated with the manufacture of surgical instruments and medical devices. Highly qualified employees adapt each individual instrument to its intended function. Therefore, most instruments are really “hand-made”.

The instruments’ finish has undergone changes lately as a result of more powerful operating lights and video camera systems being used in the OR. In spite of their superior corrosion resistance, mirror-finish instruments are available today only by special order because their light-reflecting surfaces tend to produce an irritating effect during operations. This is why state-of-the-art instruments feature a matte surface today. Whether such non-glare finish is achieved with glass beads or plastic brushes depends on the particular function of the instrument.

Final Inspection and Marking

On its way through the manufacturing cycle, each instrument undergoes a number of specified checks and tests, each of which is duly documented. Nonetheless, KLS Martin also insists on carrying out a classic final inspection. For each instrument, written inspection instructions and an inspection drawing have been put in place, against which every manufacturing batch is carefully checked. This final inspection is also documented in detail. Only those instruments that have passed this final examination are allowed to enter the next stage, where they are marked, cleaned and released for packaging and storage.

By affixing the CE-mark to the product, the manufacturer provides quality assurance. All products carrying this label fully comply with the “essential requirements” specified by the EC Directive concerning medical devices (MDD) as well as with KLS Martin’s documented and validated internal product profiles. As KLS Martin is convinced of the top quality of its surgical instruments, all instruments are covered by a lifetime warranty.

Current Reference Standards for KLS Martin Products

Standard	Content of the standard
DIN EN ISO 9001	Quality management systems - requirements
DIN EN ISO 13485	Medical devices – quality management systems - requirements for regulatory purposes
DIN EN ISO 14971	Medical devices – application of risk management to medical devices
DIN EN ISO 7153-1	Surgical instruments – metallic materials - part 1: Stainless steel
DIN EN ISO 17665-1	Sterilization of health care products – moist heat - part 1: Requirements for the development, validation and routine control of the sterilization process for medical devices
DIN 50103-3	Testing of metallic materials – Rockwell hardness testing – part 3: Modified Rockwell scales Bm and Fm for thin sheet steel
DIN 58298	Medical instruments – materials, finishing and testing
DIN 58299	Serrations for surgical instruments; profile angles, groove distances
DIN 58300	Joints for surgical instruments
DIN EN ISO/IEC 17050-1	Conformity assessment – supplier’s declaration of conformity – part 1: General requirements
DIN EN 60601-1	Medical electrical equipment – part 1: General requirements for basic safety and essential performance
MDD 93/42 Annex I	Basic requirements, annex I
MDD 93/42 Annex II	MDD 93/42 Annex II contains "EC declaration of conformity" (complete quality assurance system)
DIN EN ISO 17664	Sterilization of medical devices – information to be provided by the manufacturer for the processing of resterilizable medical devices. This information is available on the internet at http://kls-martin.com/phnet/
ASTM A967-5	Standard specification for chemical passivation treatment of stainless steel

This information is available on the internet at <http://kls-martin.com/phnet/>

General Information about the entire Instrument Range

The instruments must be cleaned, disinfected and sterilized in accordance with our preparation instructions before the first use, before each additional use and before they are returned for repair, maintenance or service. The instruments can lead to premature wear and/or risks for the patient and the user.

Ensure that the following instructions are understood and taken into account:

- Every user must read the operating instructions carefully and note them.
- All of the warnings, precautions and hazard alerts must be taken into account in particular.
- The operating instructions must be available to the user at all times. This text relates to both men and women, however references to both he and she have been omitted in order to make the text more readable.

Product Liability and Warranty

Intended Use

The instruments may only be used for their intended purpose in the specified medical fields, with use being restricted to adequately trained and qualified personnel. The treating physician or user is responsible for selecting the right instrument(s) for the surgical task / application at hand as well as for their safe handling. This includes ensuring an adequate level of training, knowledge and experience.


Warranty

Risk of damage to the instruments due to improper use!

The responsibility for proper instrument cleaning, disinfection and sterilization rests with the operator / product user. Be sure to observe your national / local regulations, including potential restrictions.

Gebrüder Martin, as the manufacturer of the products, accepts no liability for direct or consequential damage caused by improper use, handling, processing, sterilization or maintenance. Unauthorized instrument repair (by firms or persons not specifically authorized by Gebrüder Martin to perform such work) shall void the warranty given!

Non-observance of these notices, as well as improper handling or use of products supplied by us, will void your rights under the warranty. Consequently, Gebrüder Martin shall not be liable for any resulting damage in such cases.



Hotline

Should you have any questions on how to handle the unit / product or question on its clinical application, please do not hesitate to contact the Product Management:

Tel: +49 7461 706-234

Fax: +49 7461 706-312

Basic Principles of the Maintenance, Care and Preparation of Surgical Instruments

Danger of infection due to non-sterile handling!

Improper sterilization and non-sterile handling of the instruments can lead to serious health hazards for patients. All instruments must be cleaned and sterilized before using them or the first time, as well as prior to each subsequent use.

In addition, we recommend sending the products through the cleaning and disinfecting cycle three times prior to first use. Current studies have shown that such pre-treatment supports the natural formation of a passivation layer right from the start.

Processing, Cleaning, Care, Disinfection and Sterilization

Inspection and Testing Prior to Reuse

Before each use, the instruments must be thoroughly inspected for damage such as fractures, cracks or deformation, as well as for functional reliability. Special attention must be paid to cutting edges, tips, joints, box locks, ratchets and all movable parts. If wear, corrosion, deformation, porosity or other damage is detected, the instrument must be immediately withdrawn from service.

Due to their alloy, stainless steel instruments typically develop a passive film in the form of a protective layer. However, this film does not protect them well against chemical attack by chloride ions and aggressive media and liquids!

Therefore, in addition to the instrument manufacturer's endeavors to select the right materials and process them carefully, the user must make an important contribution by ensuring proper instrument processing along with adequate and regular care.

Machine Cleaning

General Notices

It is always preferable to use machines (washer-disinfectors) for instrument cleaning and disinfection because, unlike manual procedures, machine processes can be easily standardized.

Be sure to observe and follow the operating and loading instructions provided by the machine manufacturer. In addition, only the cleaning agents recommended by the manufacturer should be used for the application at hand.

- Jointed instruments should always be processed in open condition. Be sure to arrange the items so that the water can easily flow out of cannulations, blind (non-through) holes and hollow bodies.
- Complex instruments must be totally taken apart before cleaning.
- For instruments with long or narrow lumens, standard procedures should be used only if the hot disinfectant can easily flow through the lumens and safe rinsing is guaranteed.



- The instrument trays used for cleaning must always be loaded correctly to ensure proper cleaning.
- Store the instruments correctly in the tray. Be sure to prevent "rinsing shadows"!
- In accordance with DIN EN 868-8:2009 and DIN 58953-9:2010 and for ergonomic reasons, the following maximum loads are recommended:
 - Size 60 x 30 cm: 10 kg
 - Size 47 x 30 cm: 7 kg
 - Size 30 x 30 cm: 5 kg
- When removing the instruments from the machine after cleaning, be sure to check them for cleanness (visible dirt). This especially applies to cannulated instruments or those with blind holes. If necessary, repeat the cleaning cycle or pre-clean manually.
- The final rinse must guarantee that any residues left after the cleaning stage are completely removed. Be sure to use only fully demineralized water for the final rinse!

Wet and Dry Disposal

“Wet disposal” and “dry disposal” are basically different cleaning procedures.

- **Dry disposal:**
Following use in the OR, the instruments are deposited for return to the CSSD in dry condition, without applying a disinfectant or immersing them in a disinfecting solution. So for dry disposal, the instruments are not to be cleaned with or placed in a physiological saline solution.
If pre-cleaning is necessary after use, the instruments must be rinsed or wiped clean immediately.
We recommend using the “dry disposal” method for machine cleaning.
- **Wet disposal:**
Immediately after use, the instruments are to be placed into a non-fixating cleaning and disinfecting solution.
Please observe the instructions provided by the manufacturer of your disinfectant.

Detergents and Disinfectants & Temperature Ranges

- Be sure to use commercially available detergents and disinfectants that have been approved for use in each case.
- Always observe the recommended concentrations, exposure times and temperatures.
- Make sure that no residues are left on the instruments after cleaning.
- Use demineralized water for the final rinse.

	Stainless Steel Instruments	Aluminum Instruments
Cleaning/disinfection	acid / neutral / alkaline with / without the addition of tensides	neutral / mildly alkaline with / without the addition of tensides
Allowable temperature for chemical disinfection	max. 60°C / 140°F	max. 60°C / 140°F
Allowable temperature for thermal methods	max. 93°C / 199°F	max. 93°C / 199°F
Drying	at max. 125°C / 257°F	

Use of neutral-pH treating agents recommended!


Alkaline detergents and/or acid neutralizers are capable of causing discoloration – usually grayish brown – on metal-coated (e. g. TiNi) surfaces even after a short time. According to the current state of research, this does not compromise the proper functioning of the instruments. Nonetheless, using neutral-pH agents is the best option to ensure gentle treatment of your instrument surfaces.

Media containing chlorides attack surfaces!

Water with a chloride content exceeding 120 mg/l can attack the surfaces of your instruments during cleaning. If you use water with a chloride content of less than 120 mg/l, you still need to take the concentration effect during drying into account as well.

Ensure reliable drying!

Reliable drying is an essential factor for successful sterilization!



Manual Cleaning

Instruments should be disinfected and cleaned immediately after use, following the instructions given under “Cleaning and Sterilization of Surgical Instruments”. Be sure to prevent residues from drying on the instruments as this would make proper cleaning and disinfection more difficult than necessary.

The following should be observed in manual cleaning processes:

- The solutions used for manual cleaning must always be prepared in accordance with the manufacturer's instructions.
- Use a suitable brush for cleaning lumens, cannulations, blind holes and cavities, making sure that every part of the inner surface can be properly accessed.
- Use a soft brush and a neutral or mildly alkaline detergent for removing blood and other residues.
- Never use metal brushes or metal sponges for manual cleaning.
- To ensure proper instrument functioning, verify that all movable parts have been thoroughly cleaned.
- Clean jointed instruments in closed as well as open condition.
- Take instruments fully apart where applicable.
- Pay special attention to slots, ratchets, joints and box locks, narrow lumens, blind holes and other areas that are hard to access.
- Suitable trays or baskets (e. g. sterilization trays or wire baskets) must be used for storing surgical instruments properly during the cleaning process.

Ultrasonic Treatment

Effective ultrasonic cleaning requires placing the surgical instruments into suitable sterilization trays or wire baskets in an open condition. As the use of plain warm water alone cannot produce satisfactory results, it is necessary to add a suitable cleaning agent. Be sure to observe and follow the manufacturer's instructions with regard to concentration as well as temperature. Since an excessive dirt content of the cleaning solution has an adverse effect on the cleaning results, the solution must be replaced at regular intervals in accordance with the manufacturer's instructions. In the same manner, prescribed immersion or ultrasonic treatment times must be strictly observed.

As a rule, ultrasonic cleaning must always be followed by a rinsing cycle. Be sure to check the instruments (where applicable) for loosened components after the ultrasonic bath. To prevent water spots (“spotting”), fully demineralized or distilled water should be used for the final rinse.

Chemical Disinfection

- The solutions employed for chemical disinfection must always be used in accordance with the manufacturer's instructions. Never use hypochlorite-based disinfectants containing active chlorine (e. g. chlorine bleaching agent (sodium hypochlorite) or Javel water) or disinfectants containing iodine.
- Pure water must be used for preparing the dilutions specified for the chemical disinfectants. The addition of detergents is not permitted. The manufacturer's instructions regarding exposure times and concentration must be duly observed in each case.
- Disinfecting solutions must be prepared afresh every day. Extended / multiple use can easily lead to the following problems:
 - Increased concentration due to evaporation (corrosion risk)
 - Excessive dirt load (corrosion risk plus lower effectiveness).
- Following disinfection, it is important to rinse all items sufficiently under clear running water. To prevent the formation of water spots, the use of fully demineralized water is highly recommended.
- Immediately after completion of the cleaning and rinsing cycles, surgical instruments must be sufficiently dried.

Inspection after Cleaning

- Following cleaning, the instruments must be macroscopically clean, i.e. free from visible dirt or deposits.
- Instruments with stains or spots must be withdrawn from service at once and given special treatment.
- All movable parts, working tips and (scissor) blades should be inspected with particular care.
- If damage or malfunction is detected, the instrument must also be withdrawn from service immediately.

Instrument Care

“Care” means treating the instruments with instrument oil or milk (white oil-in-water emulsion). Instruments with joints or box locks (scissors, forceps, clamps, etc.) or with metal sliding surfaces (rib shears, punches, etc.) must be treated with steam-sterilizable care agents based on paraffin oil. The paraffin oil must comply with the pharmacopoeia in force at the time and must be physiologically safe as specified in the German Pharmacopoeia (“Deutsches Arzneibuch”, 10th edition (DAB 10)), “European Pharmacopoeia (Ph. Eur.)” or “United States Pharmacopoeia (USP)”.

Care agents prevent metal-on-metal friction, thus ensuring the easy movement of your instruments. Laser-marked products can be adversely affected when using basic cleaning agents containing phosphoric or hydrofluoric acid because the marking may fade away and the coding function get impaired or lost as a result.

As a rule, surgical instruments must be subjected to regular care, which means each time before a functional test is carried out. At the same time, it is important to prevent “gumming” of the joints due to an accumulative effect, especially in instruments that are continuously in use.



Sterilization

Prior to Sterilization

- Prior to sterilization, the instruments must be adequately packaged, e. g. using containers that meet EN 868-8 requirements.
- The packaging method used must comply with the relevant standards.
- Check the instruments for cleanness and integrity.
- Clean and disinfect the instruments and rinse them with distilled water, then dry them carefully,

Steam Sterilization

Danger of infection due to non-sterile handling!

Improper sterilization and non-sterile handling of the instruments can lead to serious health hazards for patients. Sterilization must be carried out according to a validated steam sterilization process, for example in a sterilizer satisfying EN 285:2009, ANSI / AAMI ST 79 requirements and validated in accordance with ISO 17665-1:2006 requirements.

It is essential to keep the steam free from any foreign bodies such as rust particles and other impurities. This helps to prevent instrument corrosion or surface damage caused by deposits. The steam used for sterilization must comply with EN 285:2009. The user instructions provided by the steam sterilizer manufacturer must be duly observed.

Instruments incorporating locks or ratchets must be sterilized in an open condition or with the ratchet set to the first notch.

The following program variants can be used:

- Validated steam sterilization process, “134°C (273°F) / 2 bar” program
- Validated steam sterilization process, “121°C (250°F) / 1 bar” program

The sterilization and holding times are subject to national provisions and regulations and therefore cannot be defined in general. It is the operator's responsibility to ensure that the desired results will be achieved with the implemented processing and sterilization procedure, including the equipment and materials used and the staff employed in the processing and sterilization department (CSSD). This requires validation and routine monitoring of the process used.

In specific, the following sources should be observed for cleaning, instrument care, disinfection and sterilization:

- ISO 17664: “Sterilization of medical devices/Information to be provided by the manufacturer for the processing of resterilizable medical devices”
- EN 285: “Sterilization - Steam sterilizers - Large sterilizers”
- ISO 17665-1: Sterilization of health care products - Moist heat - Requirements for the development, validation and routine control of a sterilization process for medical devices
- EN 556-1: Sterilization of medical devices - Requirements for medical devices to be designated “STERILE” - Part 1: Requirements for terminally sterilized medical devices
- Recommendation of the Instrument Preparation Working Group (Arbeitskreis Instrumenten-Aufbereitung / AKI), <http://www.a-k-i.org>, “Proper Maintenance of Instruments”
- Recommendation of the German Society for Hospital Hygiene (Deutsche Gesellschaft für Krankenhaushygiene / DGKH), <http://www.dgkh.de>, “Hygiene Requirements for Reprocessing Medical Devices”
- Recommendation of the Commission for Hospital Hygiene and Infection Prevention at the Robert Koch Institute (RKI) and the Springal Institute for Drugs and Medical Devices (BfArM), <http://www.rki.de>, “Hygiene Requirements for Reprocessing Medical Devices”

Information Concerning Alternative Sterilization Methods

Steam sterilization has established itself as a very safe and reliable sterilization method all over the world and therefore is considered the method of choice for goods to be sterilized that are insensitive to temperature and humidity. To this end, reference is usually made to steam sterilization using a validated steam sterilization process (see also EN 554). Consequently, there is no necessity to use other

sterilization methods – e. g. low-temperature plasma (LTP) sterilization or formaldehyde or ethylene oxide – for sterilizing steam-sterilizable medical devices. Since plasma sterilization is the subject of controversial discussion among experts with regard to its effectiveness in cavities and lumina, Gebrüder Martin does not perform validation of sterilization processes using gas/plasma sterilizers for steam-sterilizable medical devices. However, sterilization plant operators are free to validate their own alternative sterilization procedures for the medical devices to be sterilized.

Storage and Transportation

- Store the instruments in a clean, cool and dry place.
- Protect them against mechanical damage.
- Use adequate containers / packaging for safe storage and transportation.
- Handle with utmost care; never throw these products or allow them to fall down.
- Use approved sterilization packaging (complying with EN 868 / ISO 11607 requirements, for example) for sterilization and subsequent transportation and storage.
- When returning products, be sure to clean and disinfect all items and use sterile packaging.

Processing Restrictions, Disposal

Frequent reprocessing has only a minor influence on the service life of surgical instruments, which is mainly a function of use-related wear and tear or damage. Please dispose of your instruments in accordance with relevant local regulations, or have them properly recycled, once they have reached the end of their life cycle.

The national regulations for waste disposal must be complied with!

Cutting



Scissors



Bone Rongeur Forceps



Wire Cutting Pliers



Osteotomes, Chisels,
Gouges



Punches



Scalpels



Scissors

Intended use

Scissors with and without carbide inserts

Scissors by KLS Martin can be used to sever tissue, organs, bones, bandaging materials and suture materials.

Design Features



Tissue, wound edges, sutures




Surgical scissors, iris scissors




Vessel and dissecting scissors

Before using the scissors, the following points must be taken into consideration:

- The area marked must be oiled
- Visual control for impurities or surface changes
- The end parts of the scissors must be checked for breakages
- The cutting edges must be checked for integrity

Various different designs	
highMed	1 sheet standard – 1 sheet microtome blade
TC GOLD sutureCut	1 sheet standard – 1 sheet serrated blade + TC
TC GOLD	Standard blade + TC
 AQUILA TC GOLD	1 sheet standard serrated blade – 1 sheet microtome blade + TC
TCC BlackLine	1 sheet standard serrated blade – 1 sheet microtome blade + TC titanium nitride coating



One cutting edge toothed, one cutting edge with a microtome blade. Serrated blades prevent the tissue from moving, knife blades ensure a precise, atraumatic cut without the tissue being squeezed.


Wire Cutting Pliers

Intended use

In medicine, wire cutting pliers are frequently used to cut bone wires. The maximum wire diameter to be cut can be determined on the basis of the features of the wire.

For soft wires up to \varnothing 2.8 mm
 For hard wires up to \varnothing 2.2 mm

Design Features




Replaceable carbide plates



Non-replaceable carbide plates

Before using the wire cutting pliers, the following points must be taken into consideration:

- The areas marked must be oiled
- The cutting blades must be checked for damage
- If necessary replace carbide cutting blades
- A test cut can be carried out with Kirschner wire



Punches

Our range includes various different punches, including punches which cut upwards and downwards.



Cut upwards



Cut downwards



St Stainless steel



SB Solid Black
Hard-coated surface
(for extended life)





Gold spring = thin base plate

Intended use

KLS Martin punches are primarily used in neurosurgery and orthopedics for the preparation of cartilage and bones in the spine. However, the swages can also be used in other areas of medicine.

Design Features



The punches must be disassembled before cleaning using a simple, rapid disassembly mechanism. This ensures 100% cleaning of the punches.

Before using the punches, the following points must be taken into consideration:

- The area marked must be oiled
- Visual control for impurities or surface changes
- The end parts of the punches must be checked for breakages
- The punch must open and close easily
- The cutting edges must be checked for integrity
- A test cut must be carried out using a special KLS Martin silicone strip



In order to ensure 100% cleaning of our punches, special cleaning and storage trays are available for our demountable punches.

Bone Rongeurs Forceps


Intended use

Bone rongeur forceps are used to prepare parts of the cartilage and bone. They are also used to sever fine bones. They are primarily used in orthopedics.

Design Features



Standard design



Bone rongeur forceps with double gearing

Before using the bone rongeurs/forceps, the following points must be taken into consideration:



- The areas marked must be oiled
- Visual control for impurities or surface changes
- The end joints must be checked for breakages
- Checking whether all of the necessary screws are still present
- The cutting edges must be checked for integrity
- Checking whether the working ends close in parallel
- A test cut must be carried out using a special KLS Martin box (the front 2/3 of the working end must cut cleanly)

Osteotomes Chisels Gouges


Intended use

Chisels, osteotomes and gouges are used to prepare bones.


Design Features



Working part chisel



Working part osteotome



Working part gouge

Before using the chisels, osteotomes and gouges the following points must be taken into consideration:

- The cutting edges may not have any nicks
- Visual control for impurities or surface changes

Scalpels

Intended use

Scalpels are used to sever tissue, vessels and organs. A sharp, even grinding of the blade is important for good cutting ability. The grinding is either hollow or flat.

Design Features



Before using the scalpel, the following points must be taken into consideration:

- The cutting blades must be checked for damage
- Visual control for impurities or surface changes



In order to dispose of the single-use scalpel safely, please use KLING-EX (10-199-00-01) specially developed by KLS Martin for this purpose. This can be used to dispose of the blades quickly and safely.



Holding/Gripping



Clamps



Forceps



Needle Holders



Pliers

Clamps


Intended use

Clamps are instruments used to grip, hold and clamp blood vessels, tissue, organs and medical supplies.

They are classified as follows:


- Hard gripping clamps
- Soft gripping clamps
- Atraumatic clamps
- Anatomical clamps

Design Features





Atraumatic clamps:




 De Bakey Toothed



 De Bakey baby Toothed



 Cooley Toothed

Before using the clamps, the following points must be taken into consideration:

- The area marked must be oiled
- Visual control for impurities and surface changes
- The toothings must be checked for integrity
- The end joints must be checked for breakages
- The clamp must close correctly on the working part
- The lock must lock in place and may not unlock independently




Needle Holders



Intended use

Needle holders are used to grip and hold surgical suture needles during an operation. In medicine, sutures are only inserted using a needle holder. The right choice of needle holder is important. The appropriate needle holder can be determined using the list on the left hand side.

Design Features



The following carbide inserts are used for **TC GOLD** needle holders

	Division	recommended for suture material
	normal 0.5 mm	 0.6 – 4.0
	mini 0.4 mm	 4.0 – 6.0
	micro 0.3 mm	 5.0 – 12.0
	flat --	 5.0 – 12.0

Before using the needle holder, the following points must be taken into consideration:

- The area marked must be oiled
- Visual control for impurities and surface changes
- The end joints must be checked for breakages
- The carbide inserts must be checked for integrity
- The lock must lock in place and may not unlock independently
- The carbide inserts must be able to be closed into one another

Forceps



Anatomical



Atraumatic
(e.g. De Bakey toothings)



Surgical
(mouse teeth)


Intended use

Forceps are used to hold and grip tissue, organs, medical supplies and materials.

Forceps are divided into three main categories:

- Anatomical forceps
- Atraumatic forceps
- Surgical forceps

Design Features



Before using the forceps, the following points must be taken into consideration:

- Visual control for impurities and surface changes
- Check working ends for integrity
- In the closed position, the working ends must lie on top of one another perfectly
- The spring part must be checked for fissures

Pliers


Intended use

Flat pliers are used to remove bone wires.

Flat pliers are available in two different designs:

- Standard design
- Flat pliers with carbide inserts

Design Features



Standard design



Flat pliers with carbide inserts

TC GOLD

Before using the flat pliers, the following points must be taken into consideration:

- The areas marked must be oiled
- Visual control for impurities or surface changes
- The end joints must be checked for breakages
- The carbide inserts must be checked for integrity

If there is excessive wear, the carbide inserts can be replaced by an expert from the KLS Martin team.



Suctioning/Rinsing



Suction Cannulas



Cannulas



Trocars



Syringes


Suction Cannulas

Intended use

Suction cannulas are used in medical technology to suction various different liquids, such as

- Secretion
- Blood
- Fat
- Saline solution from the operating area.

Design Features



Non-demountable instruments with rinse connection must be sufficiently rinsed with detergent/disinfectant cleaning solution. Sufficient throughflow must be ensured!

Before using the suction cannula, the following points must be taken into consideration:


- The tip must be removed and the holes examined for blockages
- Visual control for impurities or surface changes
- All connections must be checked for fissures or leaks
- The shaft must be checked for deformities

Trocars

Intended use

In minimally invasive surgery, a sharp or blunt trocar is used to create access to a body cavity, and this is kept open with a tube. Various different instruments can be inserted into the body cavity through this tube.

Design Features



Before using the trocars, the following points must be taken into consideration:

- The tip must be checked for damage
- Visual control for impurities or surface changes

Cannulas

Intended use

A cannula is a hollow needle which is used to inject or tap patients with fluids or medication using a syringe. The end of a cannula is mostly filed at an angle in order to make a small cut when entering the tissue.

Design Features



Luer Lock connection cannula



LL connection syringe

Cannulas are used with a Luer Lock connection as standard. This makes it easy to attach them to a syringe.

The very sharp tip of the cannula means there is a high risk of injury.

Before using the cannulas, the following points must be taken into consideration:

- The tip must be checked for damage
- Visual control for impurities or surface changes

Syringes

Intended use

A syringe is made up of a cylindrical cavity, a piston which can move in said cavity and a working part. There are two different designs for the working part, one cone-shaped nozzle and one version with a screw thread (Luer Lock). The working part is used to administer (inject) liquid medication.

Design Features



Various different adapters for syringes



Tube connector



Catheter attachment



One-way cock

In order to guarantee safe cleaning of syringes, the piston and the cylinder must be disassembled.

After the assembly and before the use of the syringes, the seal on the syringe must be checked.

Before using the syringes, the following points must be taken into consideration:

- The syringe must be checked for damage
- Visual control for impurities or surface changes
- The function of the Luer Lock connection must be ensured



Spreading




Wound Spreaders



Rib Retractors



Abdominal Retractors



Speculums




Wound Retractors

Wound Spreaders

Intended use

Wound spreaders are surgical instruments which are used to keep an operating area open. Unlike tissue retractors, wound spreaders are self-holding instruments. This is mostly achieved using a rest lock.

Design Features



Special design with central valve



Wound spreaders are available in the following designs:

- Blunt
- Sharp
- Semi-sharp
- Full-blade
- Rigid working parts
- Movable working parts (with a joint)

Before using the wound spreader, the following points must be taken into consideration:


- The areas marked must be oiled
- The end joints must be checked for breakages
- The function of the rest lock must be checked
- The working ends (teeth) must be checked

Abdominal Retractors

Intended use

Abdominal retractors are surgical instruments which are used to keep an operating area open. Unlike tissue retractors, abdominal blades are self-holding instruments. This is mostly achieved using a rest lock.

Design Features



Special design with a central valve



Abdominal retractors are available in the following designs:

- With rigid working parts
- With movable working parts
- Special design with a central valve

Before using the abdominal retractors, the following points must be taken into consideration:



- The areas marked must be oiled
- The end joints must be checked for breakages
- The function of the rest lock must be checked
- The leaves must be checked for burrs
- The working ends must be checked for movement

Wound Retractors

Intended use

Wound retractors are used to hold tissue, organs and bones and to spread the edges of wounds.


Design Features




sharp



blunt



semi-sharp



full-blade

Tissue retractors are available in various different designs:

- Sharp
- Blunt
- Semi-sharp
- Full-blade
- One-pronged
- multi-pronged

Before using the wound retractor, the following points must be taken into consideration:


- Visual control for impurities and surface changes
- Check working ends for integrity

Rib Retractors

Intended use

In medicine, rib retractors are used to spread the sternum during a heart operation.

Design Features



Rib retractors are available in various different materials:

St Stainless steel

Ti Titanium

Al Aluminum

Rib retractors are available in various different materials:

- Stainless steel
- Aluminum
- Titanium

Only mild alkaline cleaners may be used in the preparation of aluminum rib retractors. The rib retractor must be completely disassembled for this.

Before using the rib retractors, the following points must be taken into consideration


- The area marked must be oiled
- Visual control for impurities and surface changes
- The working ends must be checked for movement
- The functionality of the driving lever and the toothed rod must be checked

Speculums

Intended use

A speculum is a medical examination instrument which is primarily used in gynecology and ear, nose and throat medicine.

Design Features



Before using the speculum, the following points must be taken into consideration:

- The areas marked must be oiled
- Visual control for impurities and surface changes
- Blades must be checked for burrs and damages
- The ease of movement of the lever must be examined
- The end joints must be checked for breakages
- The working ends must be checked for movement

Surgical Instruments
16,000 ways to operate



KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

Martin Nederland/Marned B.V.

Huizen · The Netherlands
Tel. +31 35 523 45 38
nederland@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
italia@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Moscow · Russia
Tel. +7 499 792-76-19
russia@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Shanghai · China
Tel. +86 21 5820 6251
china@klsmartin.com

Nippon Martin K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
nippon@klsmartin.com

Gebrüder Martin GmbH & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin France SARL

Mulhouse · France
Tel. +33 3 8951 3150
france@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 505 7838
malaysia@klsmartin.com

KLS Martin UK Ltd.

London · United Kingdom
Tel. +44 1189 000 570
uk@klsmartin.com

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
P.O. Box 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com

